



ABO blood typing using AmpliTaq Gold® DNA polymerase

Materials & Reagents Needed:

10X Primer Mix
AmpliTaq Gold® DNA polymerase (Applied Biosystems, Foster City, CA)
Gold ST*R 10X Buffer (Promega, Madison, WI)

10X Primer Mix:

SNP	Primer	Primer Sequence (5'- 3')	Conc.
261	Δ G ABO261d-F	AGGAAGGATGTCCTCGTG <u>T</u> TAC	1.5 μ M
	G ABO261G-F	ttaAGGAAGGATGTCCTCGT <u>I</u> GT G	2.6 μ M
	- ABO261-R	6FAM -GTTCTGGAGCCTGAACTGCT	4.1 μ M
526	C ABO526C-F	AGCTGTCAGTGCTGGAG <u>A</u> TGC	3.3 μ M
	G ABO526G-F	ttatGCTGTCAGTGCTGGAGG <u>A</u> GG	2.7 μ M
	- ABO526-R	6FAM -TCCACGCACACCAGGTAATC	6.0 μ M
803	G ABO803G-R	CCGACCCCCCGAAG <u>T</u> ACC	3.8 μ M
	C ABO803C-R	atatCCGACCCCCCGAAG <u>A</u> T C G	4.6 μ M
	- ABO803-F	6FAM -GAGATCCTGACTCCGCTGTT	8.4 μ M

Two alleles of a SNP at the 3' end of allele-specific primers are indicated in bold.

Additional mismatches included at the third or other position from the 3' end are underlined.

The tails inserted at the 5' end of one of the two allele specific primers are written in lower case.

PCR Mixture:

PCR Component	Volume/Sample
dH ₂ O	Add to 10.0 μ L
Gold ST*R 10X Buffer	1.0 μ L
10X Primer Mix	1.0 μ L
AmpliTaq Gold® DNA polymerase	0.5 μ L
DNA Template (1ng/ μ L)	1.0 μ L
Total	10.0 μ L

Thermal Cycling:

95°C for 11 minutes, then:

95°C for 20 seconds

59°C for 30 seconds

72°C for 30 seconds

For 30 cycles, then:

60°C for 60 minutes

4°C soak



Electropherograms of Six Common ABO Genotypes

