



# Multiplex SNaPshot for Body Fluid Identification

## Multiplex PCR

### Reagents Needed:

5 X Primer Mix

AmpliTaq Gold® DNA Polymerase (Applied Biosystems, Foster City, CA)

Gold ST\*R 10 X Buffer (Promega, Madison, WI)

### 5 X Primer Mix for Multiplex PCR:

	Target ID	Sequence (5'→3')	Conc. (uM)	Amplicon size (bp)
SE1	cg17610929	TTG TTG ATA TGT TTT GAA TTA TTA AG	3.0	174
		ATA ACT TCC CTT ATC AAC ACC AAC	3.0	
SE2	cg26763284-138d	TGA TTT ATA ATT ATT AGG GAG GGA AAT AG	1.0	105
		CCT AAA ACA ACC RAT TCC CAA C	2.0	
BL1	cg06379435	TTT ATT GGG GTA TTT TTA TTG GTT AG	10.0	157
		AAA ATA CAA CTT ACT CCT AAA CAC C	10.0	
BL3	cg08792630	TGT TTT AAG AGG ATG ATA AGG AA	3.0	220
		CCA CCT CAA TCC AAA CTA ACT ACA	3.0	
VF1	cg09765089-231d	TTG GTA GTT TTT GGA TTT TGG AG	3.0	137
		AAA CRT AAA ACR ACC CRA AC	24.0	
VF2	cg26079753-7d	TTT TGT GAG TGT GAG AGA TTT TTA AGA	2.0	176
		AAA ACC TCC AAA ACA AAA CCT CTA	2.0	
SA1	cg09652652-2d	GGG GAT TYG TTT YGT TAG GT	16.0	153
		CCA TTT CCC CCT TCC TAA AA	4.0	

### PCR Mixture:

PCR Component	Vol. (ul)
dH <sub>2</sub> O	~12.4
10 X Gold ST*R Buffer	2
5 X Primer Mix	4
AmpliTaq Gold (5 U/μL)	0.6 (3 U)
Bisulfite converted DNA	1 (~4)*
<b>Total</b>	<b>20</b>

### Thermal Cycling:

95°C for 11 minutes, then:

94°C for 20 seconds

56°C for 60 seconds

72°C for 30 seconds

for **34 cycles**, then:

72°C for 7 minutes

4°C soak

\* Please be aware that you should not use more than 1/5 PCR volume of bisulfite converted DNA when using Sigma's Imprint® DNA modification kit because, in our experience, it may cause PCR failure.

## Post-PCR Reaction

### Enzyme Purification of the PCR Product

#### Reagents Needed:

PCR product	10 $\mu$ L
ExoSAP-IT® (USB, Cleveland, OH)	2 $\mu$ L

#### Thermal Cycling:

37 °C for 45 minutes
80 °C for 15 minutes

### Multiplex SNaPshot

#### Reagents Needed:

10 X SBE Primer Mix  
5 X Sequencing buffer\_BigDye Termination (Applied Biosystems, Foster City, CA)  
SNaPshot™ Kit (Applied Biosystems, Foster City, CA)

#### 10 X SBE Primer Mix:

	Target ID	Sequence (5'→3')	Conc. ( $\mu$ M)	Length (nt)
SE1	cg17610929	(T) <sub>5</sub> CCG AAA CCC TCC CCA C	4.0	21
		(T) <sub>6</sub> CCA AAA CCC TCC CCA C	4.0	22
SE2	cg26763284-138d	(T) <sub>6</sub> CGC GTA ACG ACT ATA AAA CCC TC	0.3	29
		(T) <sub>8</sub> CAC ATA ACA ACT ATA AAA CCC TC	1.0	31
BL1	cg06379435	(T) <sub>17</sub> CCR ATA AAA CCT CAA ACR TAA AAC	40.0	41
BL3	cg08792630	(T) <sub>21</sub> CCR TAA TAA CTT CTA CCT ATA AAT AAA CCC	6.0	51
VF1	cg09765089-231d	(T) <sub>34</sub> TCC CCA AAT AAC AAA CRA CRA AAA TC	9.0	60
VF2	cg26079753-7d	(T) <sub>44</sub> CRA TCA ACT ACT ATA AAA ACA CC	9.0	67
SA1	cg09652652-2d	(T) <sub>48</sub> CCA CGA ATA AAT AAC CAC GAT AAA AC	15.0	74

#### SBE Reaction Mixture:

Reaction Component	Vol. ( $\mu$ l)
dH <sub>2</sub> O	~ 5
10 X SBE Primer Mix	1
5 X Sequencing Buffer	2
SNaPshot Reaction Mix	1
Purified PCR Product	> 1
Total	10

#### Thermal Cycling:

96 °C for 10 seconds
50 °C for 5 seconds
60 °C for 30 seconds
for 25 cycles

## Post-Single Base Extension

### Enzyme (SAP or CIP) Treatment

#### Reagents Needed:

SBE reaction product	10 $\mu$ L
SAP-Recombinant (USB, Cleveland, OH)	1 $\mu$ L

#### Thermal Cycling:

37 °C for 45 minutes
80 °C for 15 minutes

### Capillary Electrophoresis

#### Materials and Reagents Needed:

Dry heating block, water bath or thermal cycler  
3130 capillaries, 33 cm x 50  $\mu$ m (Applied Biosystems, Foster City, CA)  
Performance Optimized Polymer (POP4, Applied Biosystems, Foster City, CA)  
Matrix Standard Set DS-02 (dR110, dR6G, dTAMRA™, dROX™, LIZ® Dyes)  
(Applied Biosystems, Foster City, CA)  
Run Module GS STR POP4 (1 mL) E5  
GeneScan™ 120 LIZ™ Size Standard  
Hi-Di™ Formamide (Applied Biosystems, Foster City, CA)

#### Creating Matrix:

According to the ABI PRISM®SNaPshot™ Multiplex Kit protocol

#### Reagents Needed:

GeneScan™ 120 LIZ™ Size Standard	0.2 $\mu$ L
Hi-Di™ Formamide	10 $\mu$ L
SNaPshot product	1~2 $\mu$ L

#### Thermal Cycling:


95 °C for 5 minutes
4 °C soak

#### 3130 Data Collection Software:

Verify that you have chosen GeneScan Run module E5 and the DS-02 GeneScan Matrix Set.

Run prepared samples under the following conditions: injection time of 3 sec, electrophoresis voltage of 15 kV, collection time of 8 min, EP voltage of 15 kV and heat plate temperature of 60 °C.

## Electropherograms



Representative electropherograms of body fluid identification using multiplex methylation SNaPshot. SE1, SE2, BL1, BL3, VF1, VF2, and SA1 represent cg17610929, cg26763284-138d, cg06379435, cg08792630, cg09765089-231d, cg26079753-7d, and cg09652652-2d, respectively. Because all SBE primers were designed to be in the reverse direction, a blue peak represents the nucleotide G as a methylation signal and a green peak represents the nucleotide A as an unmethylation signal.