


# Laplex-12 PCR System Protocol

## Allelic size range for Laplex-12 system



## 5X Primer Mix for Multiplex PCR

Locus	Primer Sequence (5'→3')	Conc. (μM)
D4S2408	F <b>FAM</b> -AGG TAC ATA ACA GTT CAA TAG AAA GC	1.05
	R CAA TTC ATC CAC TGA AAT GAC TG	1.05
D7S3048	F TTG TAC CTT CAA GCT TCA CAA TG	1.15
	R <b>FAM</b> -CCC GCA GTC AAA AAT CTG ATA	1.15
Penta E	F GCT ACT CTG GAG GCT GAA ACA	5.00
	R <b>FAM</b> -TCA TGA TTG ATA CAT GGA AAG AA	5.00
Y-M175	F CCC AAA TCA ACT CAA CTC CA	5.00
	R <b>HEX</b> -TCT ACT GAT ACC TTT GTT TCT GTT CA	5.00
D11S2368	F TTT TAC AAT GAG GTG CAA GAA TG	1.50
	R <b>VIC</b> -TGG ATG AGT TAG ATG GGT GG	1.50
D3S2406	F <b>VIC</b> -CCT GGT ACT GCA GTG TTT CTT G	1.60
	R CCA CAT GGA GAG GCT TTA GAA	1.60
DYS391	F <b>VIC</b> -CAA GAC ACC CCA CCA CAG AT	4.25
	R <b>GTG</b> CCA TAG AGG GAT AGG TAG G	4.25

# Laplex-12 PCR System Protocol (Continued)

## 5X Primer Mix for Multiplex PCR (Continued)

Locus		Primer Sequence (5'→3')	Conc. (μM)
D8S1132	F	AGG AAA GGT TAG TGG CTT AAT GT	1.50
	R	<b>NED</b> -TCC CTC TCT CTT TCG AGC AA	1.50
D15S659	F	<b>NED</b> -ACC CTG AAG GCA GTA ATG GTT	1.90
	R	<b>GCT</b> TCC CAA CAT AAC ATA TTG CTT	1.90
Penta D	F	<b>NED</b> -GCA TGG TGA GGC TGA AGT AG	5.00
	R	TAG GTC ATG ATT TTG TGA TAT CTA AG	5.00
D6S1043	F	<b>PET</b> -CAA GGA TGG GTG GAT CAA TAG	1.65
	R	ATT GTA TGA GCC ACT TCC CAT	1.65
SE33	F	GAG TGA AAC TCC GTC RAA AGA A	5.00
	R	<b>PET</b> -CCC CTA CCG CTA TAG TAA CTT G	5.00

## Reagents Needed

- 5X Primer Mix for Laplex-12 system
- Gold ST\*R 10 X Buffer (Promega, Madison, WI)
- AmpliTaq Gold® DNA Polymerase (Applied Biosystems, Foster City, CA)

## PCR Mixture

PCR Component	Vol. (μL)
dH <sub>2</sub> O	5.5
5X Primer mix	2.0
Gold ST*R 10X Buffer	1.0
AmpliTaq Gold (5 U/μL)	0.5
Template DNA (0.5 ng/μL)	1.0
<b>Total</b>	<b>10.0</b>

## Thermal Cycling

95°C for 11 minutes, then:
94°C for 20 seconds
59°C for 60 seconds
72°C for 45 seconds
for 29 cycles, then:
60°C for 30 minutes
4°C soak

# Laplex-12 PCR System Protocol (Continued)

## Capillary Electrophoresis

### Materials and Reagents Needed

- Dry heating block or thermal cycler
- 3130 Capillary Array, 36 cm (Applied Biosystems, Foster City, CA)
- Performance Optimized Polymer (POP4, Applied Biosystems, Foster City, CA)
- GeneScan™ 500 LIZ™ Size Standard
- Hi-Di™ Formamide (Applied Biosystems, Foster City, CA)

### Mixture for CE run

GeneScan™ 500 LIZ™ Size Standard	0.15 µL
Hi-Di™ Formamide	10.0 µL
PCR product	1.0 µL

### Denaturation

95°C for 5 minutes  
4°C soak

### 3130 Data Collection Software

Application Type	HID
Injection Voltage	3.0 kV
Injection Time	5 seconds
Run Time	22 minutes

Download a PDF copy of  
this protocol

<http://forensic.yonsei.ac.kr/protocols.html>

# Genotyping result of 2800M control DNA

