




# Y-STR Multiplex System II

**DYS385, DYS437, DYS438, DYS439, DYS448, DYS456, DYS458, DYS635, GATA H4.1**



## 10X Primer Mix:

STR Loci	Primer Sequence (5'-3')	Conc.
DYS385	F <b>HEX</b> -AGC ATG GGT GAC AGA GCT A	3.2 µM
	R <u>G</u> CC AAT TAC ATA GTC CTC CTT TC	3.0 µM
DYS437	F <b>HEX</b> -GAC TAT GGG CGT GAG TGC AT	1.6 µM
	R <u>G</u> AG ACC CTG TCA TTC ACA GAT GA	1.6 µM
DYS438	F <b>HEX</b> -CCA AAA TTA GTG GGG AAT AGT TG	2.8 µM
	R GAT CAC CCA GGG TCT GGA GTT	2.8 µM
DYS439	F <b>6FAM</b> -TCG AGT TGT TAT GGT TTT AGG TCT	1.7 µM
	R <u>G</u> TG GCT TGG AAT TCT TTT ACC C	1.7 µM
DYS448	F <b>6FAM</b> -TGG GAG AGG CAA GGA TCC AA	7.5 µM
	R <u>G</u> TC ATA TTT CTG GCC GGT CTG G	7.5 µM
DYS456	F <b>HEX</b> -GGA CCT TGT GAT AAT GTA AGA TAG	1.0 µM
	R <u>G</u> TA GAG GGA CAG AAC TAA TGG AA	1.0 µM
DYS458	F <b>HEX</b> -GCA ACA GGA ATG AAA CTC CAA T	0.9 µM
	R GTT CTG GCA TTA CAA GCA TGA G	0.9 µM
DYS635	F <b>6FAM</b> -AGT GTC TCA CTT CAA GCA CCA AGC AC	4.3 µM
	R GCA GCA AAA TTC ACA GTT GGA AAA ATG T	4.3 µM
GATAH4.1	F <b>6FAM</b> -ATG CTG AGG AGA ATT TCC AA	2.0 µM
	R <u>G</u> CT ATT CAT CCA TCT AAT CTA TCC ATT	2.0 µM

## PCR Mixture:

PCR Component	Volume/Sample
dH <sub>2</sub> O	6.0 µL
Gold STR 10X Buffer	1.6 µL
10X Primer Mix	1.0 µL
AmpliTaq Gold (5U/µL)	0.4 µL
DNA Template (1ng/µL)	1.0 µL
Total	10.0 µL

## Thermal Cycling:

95°C for 11 minutes, then:

94°C for 1 minutes  
59°C for 1 minutes  
72°C for 1 minutes  
for 30 cycles, then:

60°C for 45 minutes  
4°C soak



# Electrophoresis on the ABI PRISM® 310 Genetic Analyzer

## Y-STR Multiplex System II

### Materials & Reagents Needed:

Dry heating block, water bath, or thermal cycler  
310 capillaries, 47cm x 50 µm (Applied Biosystems, Foster City, CA)  
Performance Optimized Polymer (POP4, Applied Biosystems, Foster City, CA)  
Flourescent Amidite Matrix Standards [6FAM™, TET™, HEX™, TAMRA™, ROX™]  
(Applied Biosystems, Foster City, CA)  
Run module GS STR POP4 (1 mL) D  
GeneScan™ 400HD ROX™ Size Standard (Applied Biosystems, Foster City, CA)  
Hi-Di™ Formamide (Applied Biosystems, Foster City, CA)

### Creating Matrix:

According to the ABI PRISM® 310 Genetic Analyzer User's Manual

### Preparing the Sample:

1. Prepare a loading cocktail by combining and mixing the 0.3 µL GeneScan™ 400HD ROX™ Size Standard and 20 µL Hi-Di Formamide per sample.
2. Vortex for 10 seconds.
3. Combine 20.3 µL of the prepared loading cocktail and 1.0 µL of the PCR product.
4. Preparing the allelic ladder, combine 20 µL of the prepared loading cocktail and 1.0 µL of the allelic ladder mix. Vortex the allelic ladder mix prior to pipetting.
5. Denature the samples and ladder by heating at 95°C for 5 minutes and immediately chill on crushed ice. Denature the samples just prior to loading.
6. Assemble the tubes in the appropriate autosampler, and place the autosampler tray in the instrument.

### 310 Data Collection Software:

Prepared the samples are run using the Run module **GS STR POP4 (1 mL) D** and a described above **matrix**.

Samples are injected for 5 seconds at 15,000 V and separate at 15,000V for 24 minutes with run temperature of 60°C.

### Genotyper Software:

Y-STR2plus.gta