


Y-miniplex I System

DYS19, **DYS389**, **DYS392**, **DYS437**, **DYS439**, **DYS448**, **DYS449**, **DYS456**, **DYS458**, **DYS635**



Reagents Needed:

10X Primer Mix
 AmpliTaq Gold® DNA Polymerase (Applied Biosystems, Foster City, CA)
 Gold ST*R 10X Buffer (Promega, Madison, WI)

10X Primer Mix:

STR Loci	Primer Sequence (5'→3')	Conc.
DYS19	F 6FAM -CTA CTG AGT TTC TGT TAT AGT	5.5 μM
	R <u>G</u> AT GGC CAT GTA GTG AGG ACA	5.5 μM
DYS389	F 6FAM -CCA ACT CTC ATC TGT ATT ATC T	5.0 μM
	R <u>G</u> TT ATC CCT GAG TAG CAG AAG AAT	5.0 μM
DYS392	F 6FAM -AAA AGC CAA GAA GGA AAA CAA A	1.4 μM
	R <u>G</u> AA ACC TAC CAA TCC CAT TCC TT	1.4 μM
DYS437	F HEX -GAC TAT GGG CGT GAG TGC AT	1.5 μM
	R <u>G</u> AG ACC CTG TCA TTC ACA GAT GA	1.5 μM
DYS439	F NED -ACA TAG GTG GAG ACA GAT AGA TGA T	1.6 μM
	R GCC TGG CTT GGA ATT CTT TT	1.6 μM
DYS448	F <u>G</u> CA GAA AGG GAG ATA GAG ACA TGG	0.7 μM
	R HEX -TCA TAT TTC TGG CCG GTC TGG	0.7 μM
DYS449	F NED -CTT GCT CTT TTT CTT TTC TCT CTT	5.2 μM
	R GCA CTC TAG GTT GGA CAA CAA	5.2 μM
DYS456	F HEX -GGA CCT TGT GAT AAT GTA AGA TAG	1.2 μM
	R <u>G</u> TA GAG GGA CAG AAC TAA TGG AA	1.2 μM
DYS458	F HEX -GCA ACA GGA ATG AAA CTC CAA T	1.0 μM
	R GTT CTG GCA TTA CAA GCA TGA G	1.0 μM
DYS635	F GGC TTC TCA CTT TGC ATA GAA TC	1.4 μM
	R NED -ACC AGC CCA AAT ATC CAT CA	1.4 μM

PCR Mixture:

PCR Component	Volume/Sample
dH ₂ O	4.9 μL
Gold ST*R 10X Buffer	1.6 μL
10X Primer Mix	1.0 μL
AmpliTaq Gold (5U/μL)	0.5 μL
DNA Template (100pg/μL)	2.0 μL
Total	10.0 μL

Thermal Cycling:

95°C for 11 minutes
 96°C for 1 minute, then:

94°C for 30 seconds
 59°C for 30 seconds
 70°C for 45 seconds
 for 10 cycles, then:

90°C for 30 seconds
 59°C for 30 seconds
 70°C for 45 seconds
 for 21 cycles, then:

60°C for 45 minutes
 4°C soak

Electrophoresis on the ABI PRISM® 310 Genetic Analyzer

Y-miniplex I System

Materials & Reagents Needed:

Dry heating block, water bath, or thermal cycler
310 capillaries, 47 cm x 50 μm (Applied Biosystems, Foster City, CA)
Performance Optimized Polymer (POP-4™ Polymer, Applied Biosystems, Foster City, CA)
Fluorescent Amidite Matrix Standards [6FAM™, TET, HEX, TAMRA™, ROX™]
(Applied Biosystems, Foster City, CA)
NED™ Matrix Standard (Applied Biosystems, Foster City, CA)
Run module GS STR POP4 (1mL) D
GeneScan™ 400HD ROX™ Size Standard (Applied Biosystems, Foster City, CA)
Hi-Di™ Formamide (Applied Biosystems, Foster City, CA)

Creating Matrix:

According to the ABI PRISM® 310 Genetic Analyzer User's Manual

Preparing the Sample:

1. Prepare a loading cocktail by combining and mixing the 0.3 μL GeneScan™ 400HD ROX™ Size Standard and 20.0 μL Hi-Di™ Formamide per sample.
2. Vortex for 10 seconds.
3. Combine 20.3 μL of the prepared loading cocktail and 1.0 μL of the PCR product.
4. Preparing the allelic ladder, combine 20.0 μL of the prepared loading cocktail and 1.0 μL of the allelic ladder mix. Vortex the allelic ladder mix prior to pipetting.
5. Denature the samples and ladder by heating at 95°C for 5 minutes and immediately chill on crushed ice. Denature the samples just prior to loading.
6. Assemble the tubes in the appropriate autosampler, and place the autosampler tray in the instrument.

310 Data Collection Software:

Prepared the samples are run using the Run module **GS STR POP4 (1mL) D** and a described above **matrix**.
Samples are injected for 5 seconds at 15,000 V and separate at 15,000 V for 22 minutes with run temperature of 60°C.

Genotyper Macro:

Y-miniplex I.gta

Electropherogram of Y-miniplex I Allelic Ladder and Genotyping of Standard 9948 DNA

Y-miniplex I

