


# Y-STR Analysis of Degraded DNA using In-house Massively Parallel Sequencing Panel

## Workflow



## Multiplex PCR

### Reagents Needed

5X KplexSeq-Y24 Primer Mix  
 AmpliTaq Gold® DNA Polymerase (ThermoFisher Scientific, Waltham, MA)  
 Gold ST\*R 10X Buffer (Promega, Madison, WI)

### PCR Mixture

PCR Component	Vol. (µl)
Nuclease-Free Water	11.0
10X Gold ST*R Buffer	2.0
5X Primer Mix*	4.0
AmpliTaq Gold (5U/µl)	1.0
Degraded DNA (< 100pg/µl)	2.0
<b>Total</b>	<b>20.0</b>

\*5X primer mix can be divided into 2 ~ 3 sets

### Thermal Cycling

95°C for 11 minutes, then:

94°C for 20 seconds  
 60°C for 60 seconds  
 72°C for 45 seconds  
 for x 27 cycles, then:

72°C for 5 minutes  
 4°C soak

# Y-STR Analysis of Degraded DNA using In-house MPS Panel (continued)

## Bead-based Purification of the 1<sup>st</sup> PCR product

### Materials and Reagents Needed

Agencourt® AMPure® XP beads (Beckman Coulter, Indianapolis, IN)  
Freshly prepared 80% ethanol  
DynaMag™-2 Magnet (ThermoFisher Scientific, Waltham, MA)  
Elution buffer or Nuclease-Free Water  
Dry heating block (for optional incubation in 37°C)

### Protocol

1. Remove the AMPure® XP beads from storage and let it stand for at least 30 minutes for equilibration to room temperature
2. Vortex the AMPure® XP beads to homogenize the suspension before use
3. Add 30µl of low EDTA TE buffer or Nuclease-Free Water to each sample tube that include 20µl of the 1<sup>st</sup> PCR product, and mix by tapping
4. Add x1.5 well-mixed AMPure® XP beads to 50µl each sample, then carry out procedures according to manufacturer's instructions
5. Elute PCR product in 20µl of the elution buffer or Nuclease-Free Water

# Y-STR Analysis of Degraded DNA using In-house MPS Panel (continued)

## Indexing PCR

### Reagents Needed

Nextera® XT v2 index kit (Illumina®, Inc., San Diego, CA)  
AmpliTaq Gold® DNA Polymerase (ThermoFisher Scientific, Waltham, MA)  
Gold ST\*R 10X Buffer (Promega, Madison, WI)

### PCR Mixture

PCR Component	Vol. (µl)
Nuclease-Free Water	0.4
10X Gold ST*R Buffer	3.0
Index 1 (i7)	3.0
Index 2 (i5)	3.0
AmpliTaq Gold (5U/µl)	0.6
<b>Purified</b> 1 <sup>st</sup> PCR product	<b>20.0</b>
<b>Total</b>	<b>30.0</b>

### Thermal Cycling

95°C for 15 minutes, then:

94°C for 20 seconds

61°C for 30 seconds

72°C for 45 seconds

for x 16 cycles, then:

72°C for 5 minutes

4°C soak

## Quality Control

### Materials and Reagents Needed

Agilent 2100 Bioanalyzer (Agilent Technologies, Inc., Santa Clara, CA)  
Agilent DNA 1000 kit (Agilent Technologies, Inc., Santa Clara, CA)

### Protocol

According to the Agilent 2100 Bioanalyzer and DNA 1000 kit User's Manuals  
Finally, normalize each library to 10ng/ul and pool them in equal volume

# Y-STR Analysis of Degraded DNA using In-house MPS Panel (continued)

## Bead-based Purification of the Pooled Library

### Materials and Reagents Needed

Agencourt® AMPure® XP beads (Beckman Coulter, Indianapolis, IN)  
Freshly prepared 80% ethanol  
DynaMag™-2 Magnet (ThermoFisher Scientific, Waltham, MA)  
Elution buffer or Nuclease-Free Water  
Dry heating block (for optional incubation in 37°C)

### Protocol

1. Remove the AMPure® XP beads from storage and let it stand for at least 30 minutes for equilibration to room temperature
2. Vortex the AMPure® XP beads to homogenize the suspension before use
3. Add x1.1 well-mixed AMPure® XP beads to 50µl pooled library, then carry out procedures according to manufacturer's instructions
4. Elute library in 50µl of the elution buffer or Nuclease-Free Water

## Quantify Libraries

### Materials and Reagents Needed

KAPA Library Quantification Kits (KAPA Biosystems, Wilmington, MA)  
AB 7500 Real-Time PCR System (Thermo Fisher Scientific, Waltham, MA)

### Protocol

According to the KAPA Library Quantification Kits User's Manual