



Multiplex-PCR System for mtDNA Sequence Analysis



Midi-primer Set (10X Primer Mix):

Multiplex	Amplicon	Primer	Sequence (5'→3')	Conc. (μM)
Midiplex I	P11	F15989	CCC AAA GCT AAG ATT CTA AT	8.0
		R16255	CTT TGG AGT TGC AGT TGA TG	8.0
	P13	F16347	TCA AAT CCC TTC TCG TCC C	2.0
		R006	GTG ATC CAT CGT GAT GTC TT	2.0
	P22	F155m	TAT TTA TCG CAC CTA CGT TCA	6.0
		R389-5t	ta(tta) ₁ CTG GTT AGG CTG GTG TTA GG	6.0
Midiplex II	P12	F16159-15t	(tat) ₅ CAT AAA AAC CCA ATC CAC AT	4.0
		R16410m	GAG GAT GGT GGT CAA GGG A	4.0
	P21	F015	CAC CCT ATT AAC CAC TCA CG	6.0
		R261	GCT GTG CAG ACA TTC AAT TGT T	6.0
	P31	F403-15t	(aat) ₅ TCT TTT GGC GGT ATG CAC TTT	12.0
		R614	TTT CAG TGT ATT GCT TTG AGG A	12.0

Reagents Needed:

10X Primer Mix
 AmpliTaq Gold® DNA Polymerase (Applied Biosystems, Foster City, CA)
 Gold ST*R 10X Buffer (Promega, Madison, WI)

PCR Mixture:

PCR Component	Final Conc.
dH ₂ O	17.4 μL
Gold ST*R 10X Buffer	2.5 μL 1X
10X primer mix	2.5 μL 1X
AmpliTaq Gold (5 U/μL)	0.6 μL 3.0 U
DNA Template	2.0 μL
Total	25.0 μL

Thermal Cycling:

95°C for 11 minutes, then:

95°C for 20 seconds
 55°C for 60 seconds
 72°C for 30 seconds
 for 40~42 cycles, then:

72°C for 7 minutes
 4°C soak

Sequence Analysis of mtDNA Control Region

Enzyme Purification of the PCR Product

Reagent Needed:

PCR Product	20.0 μ L
ExoSAP-IT [®] (USB, Cleveland, OH)	5.0 μ L

Thermal Cycling:

37 °C for 45 minutes
80 °C for 15 minutes

Sequencing Reaction Using BigDye™ Kit

Reagent Needed:

ABI PRISM[®] BigDye™ Terminator Cycle Sequencing Ready Reaction Kit (Applied Biosystems)
Single Sequencing Primer
BigDye[®] Terminator v1.1 & v3.1 5X Sequencing Buffer (Applied Biosystems)

Sequencing Reaction Mixture:

Reaction Component	1/4 rxn	1/8 rxn
dH ₂ O	4.0 μ L	4.0 μ L
5X Sequencing Buffer	1.0 μ L	2.0 μ L
BigDye RR Mix	2.0 μ L	1.0 μ L
Primer (1.6 pmol/ μ L)	2.0 μ L	2.0 μ L
Purified PCR Product	1.0 μ L	1.0 μ L
Total	10.0 μ L	10.0 μ L

Thermal Cycling:

96 °C for 10 seconds
50 °C for 5 seconds
60 °C for 2 minutes

for 30 cycles