

Body fluid identification by simultaneous analysis of DNA methylation and body fluid- specific bacteria

Ajin Choi, Kyoung-Jin Shin, Woo Ick Yang, and Hwan Young Lee

Department of Forensic Medicine, Yonsei University College of Medicine, Seoul, Korea



Body fluid ID based on DNA analysis


- **DNA methylation analysis**

- Target : tissue-specific differentially methylated regions (tDMRs)
- Application : Semen identification

- **Body fluid-specific bacterial detection**

- Target : 16S rRNA gene, 16S-23S rRNA intergenic spacer region
- Application : Saliva and vaginal fluid identification

➔ Need to **add more markers** or to **integrate existing DNA-based body fluid identification methods** for discrimination of more types of body fluids in one multiplex reaction



Experimental procedures

DNA extraction & quantification

Methylation-Sensitive Restriction Enzyme treatment


1 ng of DNA, 10 U of *Hha* I, 37°C for 30 min

PCR amplification


Enzyme treated DNA, 2.0 U Gold *Taq* pol, 5% DMSO, fluorescence-labeled primers

Thermal cycling


95 °C for 11 min 94 °C for 20 sec
 59 °C for 1 min) X 28 cycles 60 °C for 60 min
 72 °C for 30 sec




Schema of multiplex PCR




After MSRE treatment



Body fluid ID using multiplex PCR



Semen identification using multiplex PCR




Int J Legal Med (2013) 127:35–43

Multiplex PCR in the present study

Bacteria profiles for each body fluid


| Body fluid | N | Number of positive samples | | | | Number of negative samples |
|-----------------|----|----------------------------|------------------|---------------------|------------------|----------------------------|
| | | <i>L.crispatus</i> | <i>L.gasseri</i> | <i>S.salivarius</i> | <i>V.atypica</i> | |
| Blood | 20 | 0 | 0 | 0 | 0 | - |
| Saliva | 20 | 1 | 0 | 18 | 12 | 2 |
| Semen | 20 | 0 | 0 | 0 | 0 | - |
| Vaginal fluid | 14 | 8 | 9 | 0 | 0 | 1 |
| Menstrual blood | 14 | 8 | 8 | 0 | 0 | 3 |

Forensic evaluation - Sensitivity test





Forensic evaluation - Stability test

- Aged samples showed almost identical results compared with freshly obtained samples



Forensic evaluation – Casework sample test



Conclusion

- The multiplex PCR system, which allows combined use of 4 tDMRs for **USP49, DACT1, L81528 and PFN3**, and 4 body fluid-specific bacteria markers for *L.crispatus*, *L.gasseri*, *S.salivarius* and *V.atypica*, could be used to discriminate blood, saliva, semen and vaginal fluid-menstrual blood.
- The newly developed multiplex method ensures forensic applicability as well as high specificity, reliability and sensitivity, thereby facilitating more efficient body fluid identification in forensic casework.



YONSEI UNIVERSITY
COLLEGE OF MEDICINE

Thank you
for your attention



YONSEI UNIVERSITY
COLLEGE OF MEDICINE